

Typical Calibration Sequence

RunControl Initiates Calibration Run:

- Starts Front End Crates
- Launches Calibration Consumer(s)

Front End Crate Runs Local Calibration Procedure

- Fires charge injectors at various settings (e.g. QIE calibrations)
- Reads ADC output
- Computes gains, offsets, RMS's for local channels
- Forwards summary “C Banks” (minibank format) to Software EVB

Software EVB

- Collates all crates' “C Banks” into one bank, send to Calib CSL


Calibration CSL *recently has given trouble – try restart, see Ace Web pages*

- Sends summary banks to all requesting consumers
- Runs in “VIP Mode” guarantees deliver of all events (unlike physics CSL)




Calibration Consumers



- Previously launched by RunControl, AC++ modules
- Receives “C Banks” from Calibration CSL
- Do high level analysis and write results into Calibration Database
- Upon completion, spawn DBana script to make plots and post entry in eLog

Calibration data flow
NB: Raw data not normally kept



Crate CPUs







Dynamic script arguments filled at run time

Script that launches the AC++ job
owned by cdfprod account

Green: Running
Cyan: Not running



CalibConsumers Proxy Viewer				
File	Proxy			
QIE_0	Start	Stop	Kill	Mode
CLCQIE_0	Start	Stop	Kill	Mode
TOFQIE_0	Start	Stop	Kill	Mode
CESCALIB_0	Start	Stop	Kill	Mode
LED_0	Start	Stop	Kill	Mode



PAUSED

HALTED

RunControlStatus

State: ACTIVE

Partition: 1

ACTIVE

No error conditions

Activate: 2005.04.18 22:19:19

TS Status: OK

RunType: Physics AAA_CURRENT

TriggerType: PHYSICS_3_02 [1,639,563]

InhibitStatus: CLEAR

19:26:21

Transition: (none)

Run: 196737 0x30081 Section: 1561

12607711 Events (1.2A)

0 Calib Events

3159972 Events (1.3A)

● Clock ● Error ● SiL1T0 ● L2T0

Data Type: Beam data

AcceleratorStatus

Store Number 4,104

Tevatron Current 9.0 1E12

B0 Luminosity 36.8 E30

B0 Live Luminosity 36.8 E30

B0 Proton Losses 427.3 Hz

Accumulator Stack 34.1 mA

Fast Bunch Integrator PNG 8,102.3 E09

Tevatron Electron Lens Current 7.9 mA

AccelEvents 2005.04.19 10:27:54 Cldr:Bgn colliding physics

B0 Low Beta Current 1,974.2 Amps

Tevatron Energy 979.7 GEV

B0 Integrated Luminosity 2,586.4 nb-1

B0 Integrated Live Lum 2,337.9 nb-1

B0 Antiproton Losses 138.6 Hz

Outside Temperature 79.5 DegF

Fast Bunch Integrator ANG 840.9 E09

Store Duration 781 MINS

Or from command line:
setup fer
proxyViewer CalibConsumers

CalibConsumers Proxy Viewer

File Proxy

Consumer Name	Status	Action
QIE_0	Start	Start
CLCQIE_0	Start	Start
BSCQIE_0	Start	Start
QIEMINIPLUG_0	Start	Start
PLUGLASER2_0	Start	Start

Child: QIE_0|0

==Proxy
ProxyName: CalibConsumers
RequestId: 14

==Child
ChildName: **QIE_0**

==ChildDefinition
StartupCommand: /cdf/onln/code/cdfprod/calibration/consumer/runQiePed.6.1.4 0 b0server04 cdfonprd 285287 calqie
RestartMode: ON_REQUEST (1)
RestartTime: -1
RestartLimit: -1
ProcessPriority: -1

==ChildStatus
Alive: true
ProcessId: 343
RestartCounter: 118
ZombieCounter: 0
StatusCode: 0
LastRestartTime: Mon Nov 09 13:29:31 CST 2009
LastTerminateTime: Mon Nov 09 04:04:05 CST 2009
UserName: cdfcalib
NodeName: b0server03.fnal.gov

Consumer: QIEPED

File Browse Create Triggers

Consumer Flags: ☐ Permanent ☐ Required ☐ Vip

Parameter	Value
Directory	/Calib/QIE
Status	16777215
ProxyName	CalibConsumers
ChildName	QIE
RestartMode	ON_REQUEST
NodeName	b0server03
UserName	cdfcalib
StartupCommand	/cdf/onln/code/cdfprod/calibration/consumer/runQiePed.6.1.4 0 b0server04 cdfonprd 285287 calqie

Resolved arguments at run time

OK

Common Problems:

- Results do not show up on eLog and calibration web pages
- SoftEVB gets “stuck” during run or at end

Common Solutions:

- Restart SoftEVB via proxy control panel, easy to do, but usually not the problem
- Restart Calibration CSL, most common recent problem
 - Remember that there are *two* CSLs
 - Restart Calibration CSL on **b0server04** under the **cdfdaq** account
- If you still have problems, have a look at the calibration consumer log files:
 - /cdf/onln/cdfdaq/consumers/log
 - Each calibration type and run has unique log file
- Check cdfdaq and cdfcalib account disk quotas (done in ProcMon)